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
## Referenze Bibliografiche

1. Chen, Nanshan, et al. "Epidemiological and clinical characteristics of 99 cases of 2019 novel coronavirus pneumonia in Wuhan, China: a descriptive study." *The Lancet* (2020)
2. Wang, D, et al. "Clinical Characteristics of 138 Hospitalized Patients with 2019 Novel Coronavirus-infected Pneumonia in Wuhan, China." *JAMA* (2020)
3. Fei Zhou et al. Clinical course and risk factors for mortality of adult inpatients with COVID-19 in Wuhan, China: a retrospective cohort study. *The Lancet* (2020)
4. Claudio Ronco, Paolo Navalesi, Jean Louis Vincent, "Coronavirus epidemic: preparing for extracorporeal organ support in intensive care" [www.thelancet.com/respiratory](http://www.thelancet.com/respiratory), published online February 6, 2020;
5. Bottari G et al. "Multimodal Therapeutic Approach of Cytokine Release Syndrome Developing in a Child Given Chimeric Antigen Receptor- Modified T Cell Infusion". *Crit Care Expl* 2020; 2:e0071;
6. Napp LC et al. "Rationale of Hemoadsorption during Extracorporeal Membrane Oxygenation Support" *Blood Purification* 2019; 48: 203-214.
7. Brouwer et al. "Hemoadsorption with CytoSorb shows a decreased observed versus expected 28-day all-cause mortality in ICU patients with septic shock: a propensity-score- weighted retrospective study" *Crit Care* 2019; 317
8. Friesecke et al. "Extracorporeal cytokine elimination as rescue therapy in refractory septic shock: a prospective single-center study", *J Artif Organs* 2017; 20(3): 252-9
9. Lees NJ et al. Combination of ECMO and cytokine adsorption therapy for severe sepsis with cardiogenic shock and ARDS due to Pantón–Valentine leukocidin—positive *Staphylococcus aureus* pneumonia and H1N1, *J Artif Organs* 2016
10. Poli EC et al. "Clindamycin clearance during Cytosorb® hemoadsorption: A case report and pharmacokinetic study", *J Artif Organs* 2019
11. David S. et al. "Effect of extracorporeal cytokine removal on vascular barrier function in a septic shock patient", *Journal of Intensive Care* (2017) 5:12
12. Kogelmann K. et al. "Use of hemoadsorption in sepsis-associated ECMO-dependent severe ARDS: A case series." *Journal of the Intensive Care Society* 2018 0 (0) 1-8.
13. Calabrò MG et al. "Blood Purification With CytoSorb in Critically Ill Patients: Single-Center Preliminary Experience." *Artificial Organs* 2018, 0(0):1–6.
14. Brescia Renal Covid Task Force : Alberici F et al., GESTIONE DEL PAZIENTE IN DIALISI E CON TRAPIANTO DI RENE IN CORSO DI INFEZIONE DA CORONAVIRUS COVID-19, published March 18th by the Italian Society of Nephrology on [https://sinality.org/wpcontent/uploads/2020/03/COVID\\_guidelines\\_1703\\_finale.pdf](https://sinality.org/wpcontent/uploads/2020/03/COVID_guidelines_1703_finale.pdf) and by ERA-EDTA : [https://www.era-edta.org/en/wp-content/uploads/2020/03/COVID\\_guidelines\\_finale\\_eng-GB.pdf](https://www.era-edta.org/en/wp-content/uploads/2020/03/COVID_guidelines_finale_eng-GB.pdf)
15. Liang T (editor-in-Chief). Handbook of COVID-19 Prevention and Treatment, The First Affiliated Hospital, Zhejiang University School of Medicine, Compiled According to Clinical Experience.
16. <https://www.aferetica.com/aggiornamento-sul-possibile-uso-di-cytosorb-in-pazienti-affetti-da-covid-19/>
17. <https://www.aferetica.com/possibile-uso-di-cytosorb-in-pazienti-affetti-da-covid-19/>